


# "UZit" A-LAP™ (Type N) Retractor

## The First 3-D Soft Tissue Laparoscopic Tubular Structure Retractor

- UZit Retractor has a special design – a unique cut on the middle of the retractor's barrier provides new retraction capabilities - Laparoscopy Tubular Organ Retractor
- The UZit soft tissue retractor is a 3-D Soft Tissue Laparoscopic Tubular Retractor which allows surgeon to scoop and protect soft tissue over blood vessels etc during surgery



UZit Retractor through 5mm EZaxess Flexible Port

### Examples for utilization of UZit special Retractor:

- **UZit (N-Type A-LAP) is a great tool with special design which allows surgeon to simplify and improve the Laparoscopic nephrectomy procedure**  
The surgeon views the kidney and surrounding tissue with a flexible video scope while utilizing the UZit Retractor to retract the kidney, UZit Retractor provides great tool for scoping the kidney over its blood vessels.
- **UZit (N-Type A-LAP) can simplify dramatically the most common Laparoscopic surgical procedure performed for treatment of GERD which is a Nissen fundoplication**  
UZit special design allows surgeons to manipulate organs despite tubular on its way, UZit soft tissue retractor is the ultimate solution for retraction in one shot during Nissen fundoplication procedure.

